

AURIBOND 97



97.7% GOLD

2% PLATINUM



AURIBOND 97 Properties 97.7% Au • 2% Pt

Brinell Hardness37
Vickers Hardness40
Yield Strength5,000 psi
Ultimate Tensile Strength21,000 psi
Elongation50%
Fusion Temperature1900°F, 1038°C
Casting Temperature2150°-2200°F, 1175°-1200°C
Density19.2 gm/cm ³
Coeff. of Thermal Expansion15.5 x 10 ⁻⁶ /°C

For your convenience we accept:



Restoration courtesy of Beverly Hills Dental Studio, Beverly Hills, CA

The Ideal Choice for Single Units

- Bio-compatible constituents
- Attractive, rich, yellow color
- Harder than most comparable alloys
- Easier to finish than most comparable alloys
- Capable of fabrication with cervical collar
- Inside of coping capable of being sandblasted to produce natural matte shade

**Order AURIBOND 97 Now!
Call 800-422-7373**

Aurident is a member of Identalloy Council and maintains the highest manufacturing standards.



AURIDENT, INC.
Dental Alloys you can depend on.

*Visit our website at
www.aurident.com*