



NEW

AURITEX-LP 

40.0% GOLD
32.0% PALLADIUM
20.5% SILVER

- An economical PFM Alloy
- Compatible with most porcelain systems
- Excellent biocompatibility



NEW

CLARITY 

49.7% PALLADIUM
40.0% SILVER

- An economical, palladium-based, PFM alloy
- Low specific gravity yields more units per ounce



LH 

20.0% GOLD
20.0% PALLADIUM
41.0% SILVER

- An economical Type III crown and bridge gold alloy
- Easy to cast & work with
- Attractive soft yellow color



YX 

2.0% GOLD
40.0% PALLADIUM
25.4% SILVER

- An economical, fine-grain, Type IV crown and bridge gold alloy
- Pleasing soft yellow color
- Extremely tarnish resistant
- Copper-free for patients with allergy concerns

AURIDENT OFFERS YOU **MORE** Ways to Save Without Settling **FOR LESS.**

**Aurident Refining Services
turns your scrap into
cold, hard cash!**

Refining Services you can trust!

- State-of-the-art assay procedures assure accurate analysis
- Payments based on gold, platinum, palladium and silver contents
- Flexibility in methods of return (coins, credit, check)

SPECIAL OFFER

\$50 OFF LMT6/12

Offer expires 12-31-2012

Your next refining lot charge



37th Anniversary
AURIDENT, INC.
Dental Alloys you can depend on.

SAVE NOW! CALL 800-422-7373
For more ways to save, visit aurident.com



Aurident is a member of Identalloy Council and maintains the highest manufacturing standards. Restorations courtesy of M.R. Smile Designs, Newport Beach, CA