


# AURIDENT, INC.

## Stabilized Zirconia Blanks



Aurident's Delta Yttria Stabilized Zirconia is designed to produce outstanding restorations with high translucency and strength at an economical price. Only the highest purity components are used for Delta Yttria Stabilized Zirconia. Nano-sized stabilized zirconia particles are used in an exacting, proprietary manufacturing process resulting in a consistent, equiaxial, and uniform microstructure suitable for milling on a variety of machines.



**Delta Yttria Stabilized Zirconia is offered in:**

**Disks**

in 10, 12, 14, 16, 18 or 20mm thicknesses

**Blocks**

of 20 x 19 x 15mm or 40 x 19 x 15mm

### HIGH TRANSLUCENT DISKS

Flexural strength of 1,200 MPa  
98mm diameter

THICKNESS	PRICE
10mm	\$ 139.00
12mm	\$ 167.00
14mm	\$ 192.00
16mm	\$ 219.00
18mm	\$ 246.00
20mm	\$ 273.00

### HIGH STRENGTH DISKS

Flexural strength of 1,400 MPa  
98mm diameter

THICKNESS	PRICE
10mm	\$ 127.00
12mm	\$ 151.00
14mm	\$ 175.00
16mm	\$ 200.00
18mm	\$ 225.00
20mm	\$ 250.00

### BLOCKS

Available in High Strength and High Translucent Forms

BLOCK SIZE	FORM	PRICE PER BOX OF 5
20 x 19 x 15	High Strength	\$ 64.75
20 x 19 x 15	Translucent	\$ 72.50
40 x 19 x 15	High Strength	\$ 109.75
40 x 19 x 15	Translucent	\$ 122.50

For more information about Aurident Delta Yttria Stabilized Zirconia visit [www.aurident.com](http://www.aurident.com) or call Aurident at 800-422-7373