

A New Concept in PFM Restorations



CALL IT **A BREAKTHROUGH.**
CALL IT **THE WAVE OF THE FUTURE.**
CALL IT **VANGUARD.**



 **VANGUARD** Properties

51.2% PALLADIUM

48.0% SILVER

BALANCE IRIIDIUM & RUTHENIUM

Vickers Hardness75
Yield Strength37,000 psi
Ultimate Tensile Strength47,000 psi
Elongation9%
Fusion Temperature2350°F, 1285°C
Casting Temperature2550°-2600°F 1400°-1425°C
Density11.5 gm/cm ³
Coeff. of Thermal Expansion15.0 x 10 ⁻⁶ /°C (20°-600°C)

VANGUARD

Porcelain Fused to Metal Alloy System

The Vanguard System (*Patent Pending*) from Aurident takes PFM restorations to a whole new level of affordability and quality. This **New Oxide-Free, Base-Metal Free Alloy** enhances porcelain shade and color esthetics. It's **Hassle-Free** and **Cost-Effective!**

The Vanguard System offers:

- ✓ Clean, oxide-free castings.
- ✓ No black lines. **Ever.**
- ✓ Excellent tissue response.
- ✓ Outstanding porcelain bond strength.
- ✓ Esthetics comparable to all ceramic restorations.
- ✓ Low-cost all precious metal alloy.
- ✓ Up to 35 units per ounce.



See for Yourself. Order the Vanguard System from Aurident today and **Start Saving Tomorrow.**

 **AURIDENT, INC.**
Dental Alloys you can depend on.

To order call **800-422-7373**
or visit **www.aurident.com**

Aurident is a member of
Identalloy Council and
maintains the highest
manufacturing standards.

