

# For Outstanding PFM Restorations, Take a Good Look at **VISION**



**YOU'LL LIKE  
WHAT YOU SEE**



**Aurident's Iridium Grain Refined**

## **VISION**

**Over two decades of proven reliability.**

- Designed for porcelain fused to metal restorations.
- Contains Iridium, a scarce noble metallic element of the Platinum group, to insure exceedingly fine grain structure.
- High yield strength at elevated temperatures makes VISION resistant to margin creep.
- Compatible with most porcelain systems.

**TO ORDER, CALL 800-422-7373**

See our complete product line at [www.aurident.com](http://www.aurident.com)

**80% PALLADIUM  
2% GOLD  
SILVER FREE**



**Vision Properties**

2.0% Au • 0.1% Pt • 80.0% Pd

Vickers Hardness	.....265
Yield Strength	.....143,000 psi
Ultimate Tensile Strength	.....160,000 psi
Elongation	.....10%
Fusion Temperature	.....2150°F
Casting Temperature	.....2450°-2500°F
Density	.....11.1 gm/cm <sup>3</sup>
Coeff. of Thermal Expansion	...14.0 x 10 <sup>-6</sup> /°C

Restoration courtesy  
of Tooth Effects,  
Newport Beach,  
California

Aurident is a member of  
Identalloy Council and  
maintains the highest  
manufacturing standards.



**AURIDENT, INC.**  
*Dental Alloys you can depend on.*



ISO 13485:2003  
FM 91564