

WX

Quality. Savings. Aurident's WX Delivers the Best of Both Worlds.



57.5% GOLD

0.1% PLATINUM

2.5% PALLADIUM

28.1% SILVER



WX Physical and Mechanical Properties

Vickers Hardness	.150 (Q) / 230 (H)
Yield Strength	.60,000 psi (H)
Ultimate Tensile Strength	.115,000 psi (H)
Elongation	.19%
Fusion Temperature	.1725°F, 940°C
Casting Temperature	.1800°-1850°F, 980°-1010°C
Density	.13.8 gm/cm ³

Aurident's WX High Noble Gold Casting Alloy lowers cost without compromising quality. WX offers dental laboratories:

- A versatile, gold casting alloy with a **BRILLIANT YELLOW COLOR**
- Exceptional tarnish resistance
- Ease of casting and workability
- Outstanding mechanical properties
- The ideal alloy for all Type III or Type IV crown and bridge applications

For superb quality made more affordable, get the best of both worlds... WX from Aurident.

To order Call 800-422-7373.

For your convenience we accept:



AURIDENT, INC.
Dental Alloys you can depend on.

Visit our website at www.aurident.com

Aurident is a member of Dentalloy Council and maintains the highest manufacturing standards.

Restorations courtesy of Glidewell Laboratories, Newport Beach, CA