



3rd Party Material Packages available for AccuFab Printers (D1s,L4D)

Date: 9 Jun,2022

Brand	Material	Layer Thickness	AccuFab-D1s	AccuFab-L4D
Nextdent	C&B MFH N1	0.05	√	√
	C&B MFH N1.5	0.05	√	√
	C&B MFH N2	0.05	√	√
	C&B MFH N2.5	0.05	√	√
	Denture 3D+ Dark pink	0.05	√	√
	Denture 3D+ Light pink	0.05	√	√
	Denture 3D+ Opaque pink	0.05	√	√
	Denture 3D+ Red pink	0.05	√	√
	MFH-N1.5	0.05	NA	√
	MFH-N2.5	0.05	NA	√
	Ortho Flex	0.1	√	√
Keystone	KeyCast	0.05	√	√
	KeyGuide	0.1	√	NA
	KeyMask	0.1	√	√
	KeyModel	0.1	√	√
	KeyModel Ultra Ivory	0.05	√	√
		0.1	NA	√
	KeyModel Ultra Light Grey	0.05	√	√
		0.1	NA	√
	KeyOrtho IBT	0.1	√	√
	KeyOrtho Model	0.1	√	√
	KeySplint Hard	0.1	√	√
KeySplint Soft	0.1	√	√	
KeyTray	0.1	√	NA	
pro3dure	GR-10 guide	0.1	√	NA
	GR-11 tray	0.1	√	NA
	GR-12 cast	0.05	√	NA
	GR-14.1 denture	0.075	√	√
	GR-17 temporary A1	0.05	√	√
	GR-17.1 temporary lt	0.05	√	NA
	GR-18.1 IB	0.1	√	√
	GR-19 OA	0.1	√	√

	GR-21 Try-In	0.05	√	NA
	GR-22 flex	0.1	√	√
HARZ Labs	DENTAL SAND A1-A2	0.05	√	√
	DENTAL BLEACH	0.05	√	√
	DENTAL PINK	0.05	√	√
DETAX	FREEPRINT cast 2.0 385	0.05	√	√
	FREEPRINT denture 385 translucent pink	0.05	√	√
	FREEPRINT gingiva 385	0.05	√	√
	FREEPRINT IBT 385	0.1	√	√
	FREEPRINT model T 385	0.1	√	√
	FREEPRINT ortho 385	0.075	√	√
	FREEPRINT ortho 405	0.05	NA	√
		0.075	√	NA
	FREEPRINT splint 2.0 385	0.1	√	√
	FREEPRINT temp 385 A2	0.05	√	√
FREEPRINT tray 2.0 385	0.1	√	√	
dentona	optiprint clara transparent	0.1	√	√
	optiprint laviva light pink	0.1	√	√
	optiprint lumina A2	0.05	√	√
	optiprint zero red	0.05	√	√
Liqcreate	Dental Model Pro Beige	0.05	NA	√
		0.1		√
	Dental Model Pro Grey	0.05		√
		0.1		√
	Gingiva Mask	0.05		√
		0.1		√
	Premium Model	0.05		√
		0.1		√
	Wax Castable	0.025		√
		0.05		√
0.1		√		
Gorky Liquid	Dental Castable	0.05	√	√
	Dental Crown A1A2	0.05	√	√
	Dental Model	0.05	√	√

	Dental Flex	0.1	√	NA
Dreve	FotoDent IBT	0.1	√	√
	FotoDent tray	0.1	√	√
BASF	Ultracur 3D EL150	0.05	√	NA
	Ultracur 3D FL60	0.075	√	
	Ultracur 3D RG35	0.075	√	
	Ultracur 3D RG50	0.05	√	
		0.05	√	
	Ultracur 3D ST45B	0.075	√	
		0.1	√	
	Ultracur 3D ST45M	0.05	√	
		0.075	√	
		0.1	√	
	Ultracur 3D ST80	0.05	√	

Related Articles

- **Shining3D Dental Material Certificates**
(<https://support.shining3ddental.com/en/support/solutions/articles/60000707063-shining3d-dental-material-certificates>)
- **Post Processing for Shining3D Material**
(<https://support.shining3ddental.com/en/support/solutions/articles/60000700409-post-processing-for-shining3d-material>)
- **Why the material color looks different?**
(<https://support.shining3ddental.com/en/support/solutions/articles/60000701672-why-the-material-color-looks-different->)